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# 1. IMPORTANT SAFETY INSTRUCTIONS

- The thermometer is suitable exclusively for measurement of the body temperature.
- Always respect the duration of measurement until the signaling tone [beep] must be maintained without exception!
- The thermometer contains small parts (battery, etc.) that can be swallowed by children. Keep the thermometer out of the reach of children.
- Damage may occur if the thermometer is dropped or by striking on a hard surface.
- Avoid ambient temperatures above 60°C. NEVER boil the thermometer!
- Use only the commercial disinfectants listed in the section «cleaning and disinfecting» to clean the thermometer. The thermometer must be intact when immersed in liquid disinfectant.
- Warning:** The thermometer contains a long-life battery. For disposal, please comply with local regulations for special waste.
- The manufacturer recommends that the accuracy is verified by an authorized laboratory every 2 years. Please contact our Customer Help line (0845 1308050) for advice (for UK only).

## 2. SWITCHING ON

Press the on/off button ❸. A short beep-signal sounds, the device is on. All symbols show on the LCD display ❷ for a few seconds as a means of checking that the thermometer is functioning correctly. If available the last correctly measured temperature is displayed indicated by the "M" symbol "by pushing and holding the on/off switch for more than 3 seconds. After this the device is ready for measurement indicated by a flashing "°C" symbol and a flashing "L" at an ambient temperature of less than 32°C.

Proper functioning of the thermometer is tested automatically each time it is turned on. If a malfunction is detected, this is indicated by «ERR» on the display, and a measurement becomes impossible. In this case, the thermometer must be replaced.

## 3. HOW TO USE THE THERMOMETER

Position the sensor tip ❶ in the preferred place of measurement – mouth, rectum or armpit.

## In the mouth (oral)

Place the thermometer under the tongue. The measuring sensor ❶ must be in good contact with the tissue. Close your mouth and breath evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air. Approx. measuring time: 1 minute. Normal temperature between 35.7 °C and 37.3 °C.

## In the anus (rectal)

This is the most reliable measuring method, and is especially appropriate for infants and small children. Carefully insert the tip of the thermometer 2 to 3 cm into the anal aperture. Approx. measuring time: 1 minute. Normal temperature between 36.2 °C and 37.7 °C.

## In the armpit (axillary)

From a medical viewpoint, this method will provide inaccurate readings, and should not be used if precise measurements are required. Approx. measuring time: 1½ minutes. Normal temperature between 35.2 °C and 36.7 °C.

If the signaling 4-tone (beep-beep-beep-beep) sounds 10 times and the «°C» is no longer flashing, this means that the measurement is finished and that the thermometer is ready for reading.

Warning! To prolong the life of the battery, turn it off after use by pressing the operating button briefly. Otherwise the thermometer will automatically turn off after about 10 minutes.

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**Omron Healthcare (UK) Ltd**